


Effectiveness of Palm Vein Recognition Biometric Technology to Prevent Patient Misidentification Errors

Rekha Mody, M.D.

Namita Gandhi, M.D.



Pilot Parameters

- 4 weeks
- Outpatient imaging only
- Regional hospital
 - Small, controlled environment
 - Quick throughput
 - Cross section of most imaging studies
 - Studies with contrast performed
 - Intermediate volume
- ★ Patient enrollment & authentication in same encounter

Metrics

Measure	Goal
# Misidentification Errors	0
Consistency	<2%
a. % enrolled patients not identified when re-scanned	
b. % patient unable to be enrolled	<3%
Enrollment %	≥ 80%
Authentication % (Point of Care)	≥90%
Patient experience	Score ≥4/5 all questions
Caregiver experience	Score ≥4/5 all questions
IT Success Assessment / EMR Integration	Qualitative
Enrollment Time Efficiency: Palm vein recognition device vs Standard process	Observational

Pilot Summary Data

Item	%/Number s
Total Unique Patients Seen	570
Total Patient Visits (includes repeat visits)	596
Total Exams (8 Unique Modalities)*	759
Total Enrollments (unique patients)	455
Enrollment % (enrolled patients / total unique patients)	80%

**X-ray, CT Scan, MRI, Ultrasound, Mammogram, Bone Density, Nuclear Medicine, Flouroscopy*

Pilot Summary Data

Item	Percent/Number
Authentication % (Point of Care)	75% (369)
Percentage of Enrolled Patients who enrolled on 1 st Visit	99.6% (453)
% Enrolled Patient with Repeat Encounters	3.5% (16)
% Repeat patients who were re-authenticated at the front desk	88% (14)
% Repeat patients who were re-authenticated at the POC	50% (8)
 % Accurate identification of re-authenticated repeat patients (front desk & POC)	100%

Overall Metric Assessment

Measure	Goal	Actual	Rating
# Misidentification Errors	0	0	
Consistency	<2%	3%	
a. % enrolled patients not identified when re-scanned			
b. % patient unable to be enrolled	<3%	1.5%	
Enrollment %	≥ 80%	79%	
% Authentication Point of Care	≥90%	75%	
Patient experience	Score ≥4 all questions	4.1	
Caregiver experience	Score ≥4 all questions	3.5	
IT Success Assessment / EMR Integration	Qualitative Measures	See Qualitative IT Scorecard	
Enrollment Time Efficiency: Palm vein recognition device vs Standard process	Observational	>Time Palm vein recognition device vs Standard	

Key Takeaways-Quality & Safety

- Safe
- Accurate
- Consistent
- Enrollment % nearly at target
- Authentication % below target



Key Takeaways-Experience

- Excellent to very good overall patient experience scores
- Very Good to good overall caregiver experience scores
- Better scores for “extremely private” & “extremely safe” questions with Palm vein recognition device vs Standard

Key Takeaways-Demographics

- Nearly all who enrolled did on 1st visit
- Females and Males decline at equal %
- Younger age groups declined at a 2X higher % than older
- Majority of declines “not interested”
- Very few declined due to COVID/Health concerns

Key Takeaways-Operations

- Increase time for Palm vein recognition device enrollment
 - Equal or less time at point of care
 - As we enroll more → less time at front desk
- “Green light” ease of EMR integration